

Product Name :
L-Spanner

Product Code :
ASTRALAB-SUPPLYA115002



Description :

L-Spanner

Technical Specification :

High quality carbon steel.

Chrome plated finish for rust resistance.

Essential socket wrench in every service station and garage, ideal tool for car wheel nuts.

Hexagon head 21mm

Carbon steel chrome plated

With chisel end / without chisel end

Size: 8-19' D=12.5-27.5, L=220-250 H=85-95.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265