

Product Name :
Radiation and Insolation Data Logger

Product Code :
ASTRALAB-SUPPLYA118001



Description :

Radiation and Insolation Data Logger

Technical Specification :

Pyranometer is a precision solar radiation measurement instrument with a ground glass dome, and is designed for long-term outdoor exposure.

It is an affordable and effective alternative to shade-ring pyranometers, pyrheliometers and traditional sunshine recorders.

Over accuracy total or global and diffuse radiation = \pm 5% daily integrate , \pm 5% \pm 10 w.m-2 hourly averages
 \pm 8%, \pm 10 w.m-2

Resolution: =0.6w.m2

Range: 0 to >2000w.m-2

Analogue output sensitivity: = 1wm-2

Temperature range: =-40 to 70 oc

Tilt response: = negligible

Power requirement: 10 v-15V



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265