

Product Name :
Digital Permeability Tester Machine

Product Code :
ASTRALAB-SUPPLYA122002



Description :

Digital Permeability Tester Machine

Technical Specification :

Microcomputer control system, high accuracy sensor, adjustable air volume

Equipped with large digital display of parameters and integral printer

Pressure difference ranging: 50-500pa

Measuring range: 5-10000mm/s

Measuring accuracy: $\leq \pm 3\%$

Sample thickness: $\leq 5\text{mm}$

Sample test area 5, 20, 50, 100cm²

Office plate (original nozzle, 8 pcs) diameter, mm: 5, 1, 1.5, 2.5, 4.5, 8, 14.5, 26

Power supply 220v 50/60Hz, 1300w

Dimension (L*W*H), mm: 700*1000*1000

Weight, 100kg.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265