

Product Name :
Egg Incubator

Product Code :
ASTRALAB-SUPPLYA131009



Description :

Egg Incubator

Technical Specification :

Double wall standard fabrication exterior wall made out of thick steel duly finished in powder coated paint.

Inner chamber made out of highly polished thick aluminum.

Temperature range ambient to 70oC controlled by imported thermostat with accuracy +0.5oC unit provided with manual tilting device for egg trays.

To work on 220 VAC supply.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265