

Product Name :
Cutaway Model of Motors

Product Code :
ASTRALAB-SUPPLYA132003



Description :

Cutaway Model of Motors

Technical Specification :

Sectioning of an actual industrial electrical component (popular motor brands and models are chosen for industrial relevance, including selections from Leeson, Dayton, GE, Baldor, Marathon, WEG, Siemens, Leeson, and others).

Careful sectioning of full length of the device.

13-gauge, formed-steel, powder-coated mounting baseplate.

Provision for mounting on related DAC Worldwide bench, workstation, and storage products.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265