

Product Name :
Network Analyzer

Product Code :
ASTRALAB-SUPPLYA132009



Description :

Network Analyzer

Technical Specification :

Two-port vector network analyzer with a full S-parameter test set for bidirectional measurements

Wide dynamic range of typical 120dB

Measurement bandwidths from 1Hz to 500KHz

Wide 10.1" WXGA multitouch screen

System keys for setup, presets, settings

Clearly structured user interface

Undo/redo softkey for user-friendly operation

Fully integrated context-sensitive help menu

Calibration units for quick calibration

De/embedding functionality and fixture compensation



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265