

**Product Name :**  
Sensor and Transducer Training Kits

**Product Code :**  
ASTRALAB-SUPPLYA138001



**Description :**

Sensor and Transducer Training Kits Supplier in Africa and Lab Equipment South Africa, Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Ivory Coast, Cameroon, Niger, Burkina Faso, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan.

#### Technical Specification :

Sensor and Transducer Training Kits

To study of humidity sensor transducer

Humidity transducer capacitive for 90% non condensed range 0 to 100%

RH humidity chamber

3.5 digit display voltmeter thru 20 turn pot/preset.

Sockets at different points for observations

Built in IC regulated power supply

Interconnection : 2/4mm banana sockets

220V  $\pm 10\%$ , 50Hz mains operated

Enclosed in an attractive ABS plastic cabinet with cover

Standard Accessories : User's Manual with patch cords.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265