


Product Name :
Digital Encoder

Product Code :
ASTRALAB-SUPPLYA138005



Description :

Digital Encoder

Technical Specification :

Encoder & Decoder Trainer has been designed specifically to demonstrate the functionality of Encoder & Decoder Circuits. This kit is an extremely useful laboratory kit which provides facilities for hands on experience to various experiments like study of decimal to binary encoder, binary to decimal decoder etc.

This kit comes with built-in fixed DC Power Supply, set of built-in logic level generator and indicator. This is helpful in making the life of instructor's simple, at the same time while making experimentations on it.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265