

Product Name : Unit Fluid Bed Dryer, Digital	Product Code : ASTRALAB-SUPPLYA140003
--	---



Description :
Unit Fluid Bed Dryer, Digital
Technical Specification :
High Precision and effective drying technology.
Pilot scale table top model.
Easy to install, operate, clean, convenient & user friendly.

Inbuilt design with Inlet hot air and exhaust air blower.

Turbine blower assembly for air suction.

Electric Heater of Stainless Steel with fins.

20 micron filters for Inlet Air.

Sampling port on product container.

All gasket white silicon food grades.

Manual Damper for Inlet and Outlet Air duct.

Inspection glass window for process visualization.

Safe earthing arrangement.

Filter bag MOC-Satin.

Pneumatic shaking for filter bag.

Auto timer for batch processing.

Digital Temperature Controller with sensor for inlet air.

Digital AMP indicator on front panel.

DP across the product container.

DP across the filter.

Pneumatically Operated Filter bag Shaking.

Pneumatically Operated Outlet damper.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265