

**Product Name :**  
Laboratory Scale Spray Dryer

**Product Code :**  
ASTRALAB-SUPPLYA140010



**Description :**

Laboratory Scale Spray Dryer

**Technical Specification :**

Individual controller for temperature, vacuum and feed rate using control knobs.

Spray dry with aqueous as well as solvent based solutions.

Compact table top model; enables spray drying trial with even 30 ml samples.

Designed with a SS material & reinforced with flameproof motors, to withstand the effect of acid & solvent spray drying.

Simple to operate and clean, economical and gives quick results.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265