

Product Name : Thermostatic Water Bath	Product Code : ASTRALAB-SUPPLYA141003
--	---



Description :
Thermostatic Water Bath
Technical Specification :
Made of Stainless steel sheet. Top of the bath and concentric rings are also of Stainless Steel.
Heating elements have an automatic cut off device when the bath gets dry. Standard double wall construction, inner being of stainless steel

outer of M.S. duly staving paint finish and gap filled for temp. insulation with glass wool.

Top cover with 75 mm (3?) holes with concentric rings.

Temp. range 5o C above ambient to 95o C type element with accuracy of + 1o C.

Supplied complete with plug cord.

Panel provided with pilot lamps, thermostat control knob etc.

Features:

This is provided with ISI mark 1500 watts heating device.

Removable concentric rings to accommodate containers of various sizes.

Controlled by capillary type thermostat from ambient to 110° C, within +1° C., two indicator lights & one on/off Rotary Switch.

The double walled chamber is filled with high-grade insulation to prevent heat loss.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265