

**Product Name :**  
Laboratory Autoclave 18L capacity

**Product Code :**  
ASTRALAB-SUPPLYA141005



**Description :**

Laboratory Autoclave 18L capacity

**Technical Specification :**

Easy operation, safe and reliable.

The sealing ring is made of silicon rubber.

Two safety valves keep double safety (for timing control type).

It has self-adjustment function (for timing control type)

It is well sealed and effectively prevents steam leakage.

Over temperature & over pressure auto-protection.

Electric heating or gas heating can be used (for Dual fuel Type).

The sterilizer is provided with a power switch, user can control the operation easily.

Immersion type electric heating pipe, thermal efficiency is high, stainless steel material tube is durable to use.

Steam-water inner circulation system: No steam discharge and the environment for sterilizing will be clean and dry.

Double groove pressure gauge, could clearly indicate, convenient observation, reliable performance.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265