

**Product Name :**  
Laboratory Centrifuge 15ml x 8 Bucket

**Product Code :**  
ASTRALAB-SUPPLYA141006



**Description :**

Laboratory Centrifuge 15ml x 8 Bucket

**Technical Specification :**

Economical unit for small capacity centrifuging requirements.

The standard unit is fitted with resiliently mounted motor with 3500 RPM for vibration free performance, a built in 5 step speed regulator and one non removable angle motor.

Complete with tube carriers.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

