

**Product Name :**  
Flow Test Training Kit

**Product Code :**  
AU - TR19



**Description :**

Flow Test Training Kit, We are manufacturer, supplier, exporter and solution provider in quality Flow Test Training Kit

**Technical Specification :**

Flow Test Training Kit - uses reconditioned engine using Mitsubishi, Mazda, Honda, Toyota or any other models that are available at the time. Some motor engines comes with a Simulator Engine Test Board, which is used for testing, measuring or simulating any signals of sensors and to activate the actuators of the engine, which is connected from the terminals of Engine Control Unit (ECU).

The training kits display use real car parts or systems and are operated either in manual or automatic.

Some of the main components displays are:

- Venturi meter
- Orifice meter

- Pilot static tube
- Flow meter with manometer for monitoring pressure



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265