

Product Name : Multi-Board (MB-85 & PLCC)	Product Code : BRD-LABAST-006
---	---



Description :
Multi-Board (MB-85 & PLCC), We are manufacturer, supplier, exporter and solution provider in quality Multi-Board (MB-85 & PLCC)
Technical Specification :
<p>Multi-Board (MB-85 & PLCC)K&Hs PLCC test boards offer users a fast, safe and solderless method of prototyping circuits with surface mount chip carriers, eliminate time-consuming soldering and expenses of building permanent connections. The PLCC test boards are composed of a chip carrier socket and special designed pins with high quality PC board. The PLCC test board comes in 5 sizes : 20, 28, 44, 68 and 84 pins.</p> <p>Combine PLCC test boards with MB-85 multiboard, a convenient and useful tool will be created. The MB-85 multiboard's center channel distance can be adjusted according to IC size. Aside from PLCC test boards, LSI, MSI and DIP size chips can all be built easily.</p> <p>MB-85</p> <ul style="list-style-type: none">- Easy and fast prototyping with surface mount chip carriers- Eliminates time-consuming soldering- Compatible with K&H solderless breadboard <p>SPECIFICATIONS</p> <ul style="list-style-type: none">- 5 interconnected clip.....128- 25 interconnected clip.....8- IC capacity (14 Pin).....9

- Weight (gm).....140

PLCC

Model

A

B

C

D

Weight

PL-20

25.4(1.0")

22.86(0.9")

31.88(1.3")

27.06(1.1")

5.6g

PL-28

27.94(1.1")

33.02(1.3")

34.42(1.4")

37.22(1.5")

7g

PL-44

33.02(1.3")

53.34(2.1")

39.5(1.6")

57.54(2.3")

12g

PL-68

40.64(1.6")

83.82(3.3")

47.12(1.9")

88.02(3.5")

23g

PL-84

45.72(1.8")

104.14(4.1")

52.2(2.1")

108.34(4.3")

30g



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265