

**Product Name :**  
Bench grinders

**Product Code :**  
AS3E25



**Description :**

Bench grinders Grinder 12" 180 Grit, 3600 rpm;  
- Two 12" Platen (120 or 180 grit,  
variable speed

The compact and Bench Grinder Machine for stationary grinding, deburring and finishing operation.  
Absolutely different with a standard bench grinder. Ideal for 200mm Out dia and width up to 50mm valgro  
deburring and finishing wheels.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265