

**Product Name :**  
Gyroscope

**Product Code :**  
TM 630



**Description :**

Gyroscope Standard: Brass wheel is 7.5cm diameter and 3cm at its rim carefully balanced and all the three pairs of pivots are adjustable

- Complete fitted on metal base.
- Bench space: 1.2m x 0.6m, speed control unit: digital tachometer and stop watch.
- Dimensions: 650mm x 510mm x 340mm, Nett weight:22Kg



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265