

Product Name :
Torsional Vibration Apparatus

Product Code :
TM 163



Description :

Torsional Vibration Apparatus Specification

- [1] Student experiment on torsional vibration
- [2] Torsion wires, made of aluminium
- [3] Mass disc made of galvanised steel
- [4] Mass ring made of galvanised steel
- [5] Strong chuck for fitting the wires
- [6] Wall mounting plate made of anodised aluminium

Technical Data

Torsion wires

- length: 1100mm

- diameter: 2mm, 3mm, 4mm, 5mm and 6mm

Mass disc

- diameter: 160mm

Mass ring

- outer diameter: 160mm

- inner diameter: 100mm

Dimensions and Weight

l x w x h: 200 x 150 x 270 mm

Weight: approx. 11 kg

Scope of Delivery



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265