

Product Name :
Vibration of Spiral Spring Apparatus

Product Code :
ASgd409



Description :

Vibration of Spiral Spring Apparatus Spiral spring

- stiffness: $c=238\text{Nmm}$
- cross-section: 10x1mm
- effective length: approx. 800mm
- inside radius: 10mm
- outside radius: 50mm

Weights made of stainless steel

- mass: 0.5kg
- diameter: 50mm

Levers made of stainless steel

- diameter: 8mm
- length: 154mm
- arm radius of the rotating mass: 36...150mm

Angle measuring scale

- 0...360°
- resolution: 1°

Base plate

- w x h: 200x250mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265