

Product Name :
Marcet Boiler

Product Code :
HE 169



Description :

Marcet Boiler SOLTEQ MODEL: Pressure Vessel: 3 litre capacity made of stainless steel, Pressure gauge: 0 to 20 bar (Bourdon type), Thermometer: Glass, 0 to 200 °C, Transducers: Temperature and pressure, Operating Pressure: Max. 15 bar, Heater: 2 kW, immersion type. DIGITAL INSTRUMENTATIONS: 2 unit of digital indicator, 1 unit of pressure transmitter, SOLDAS DATA ACQUISITION SYSTEM: A PC with latest Pentium Processor, An electronic signal conditioning system, Standalone data acquisition modules, windows software for data logging, signal analysis, display tabulation and graphing.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265