

Product Name :
TILTING SEDIMENT FLOW CHANNEL, 600 mm. wide

Product Code :
FM161



Description :

60 cm wide x 70 cm high x 18 m long, -1 to +3% motorized tilting ; can be used for study of flow phenomena without sediment or study of sediment transport and bed form in open channel ; complete with pump, feed tank, storage tank, sluice gates, motorized trolley, sand, sediment feeder, sediment trap, accessories and models.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265