

**Product Name :**  
TILTING SEDIMENT FLOW CHANNEL, 600 mm. wide

**Product Code :**  
FM161



**Description :**

60 cm wide x 70 cm high x 18 m long, -1 to +3% motorized tilting ; can be used for study of flow phenomena without sediment or study of sediment transport and bed form in open channel ; complete with pump, feed tank, storage tank, sluice gates, motorized trolley, sand, sediment feeder, sediment trap, accessories and models.



## **Astra Scientific**

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265