

**Product Name :**  
Cable Tensiometer

**Product Code :**  
AEE 4



**Description :**

Cable Tensiometer The Cable Tensiometer is calibrated for use on standard, multiple-strand, nonjacketed, aircraft cable, made of stainless steel or carbon steel with zinc or tin over zinc coating (MIL-DTL-83420 Specification, type I, comp. A and comp. B and Boeing BMS 7-265) of 7 x 7 stranding in diameters of 1/16 and 3/32 inch, and 7 x 19 stranding for 1/8, 5/32, 3/16, 7/32, and 1/4 inch.  
One analogue the other digital.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265