


**Product Name :**  
Concrete Cube Mould

**Product Code :**  
ASTConcrete-0004



**Description :**

Concrete Cube Mould

**Technical Specification :**

We are manufacturing Cube Moulds as per international standards (BS and IS). These cube moulds are available in different sizes. For the metric size cube mould, the faces are machined flat to  $\pm 0.2\text{mm}$  accuracy and finished to within  $0.2\text{mm}$ . For the inch size moulds, the faces are machined flat to  $\pm 0.01$  inch to within  $0.01\text{in}$ . All cube moulds are supplied complete with the base plate.

- Cube Mould (Cast Iron - for Mortar) 50mm, Single Gang Cube Mould (Cast Iron - for Mortar) 50mm, Three Gang Cube Mould (Mild Steel - for Cement) 70.6mm, with Loose Base Plate Cube Mould (Cast Iron - for Concrete) 75mm Single Gang Cube Mould (Cast Iron - for Concrete) 100mm Light Weight Cube Mould (Cast Iron - for Concrete) 100mm Clamp Type (8.5 KG) Four-Part Cube Mould (Cast Iron - for Concrete) 150mm Light Weight 8KG Cube Mould (Cast Iron - for Concrete) 150mm ISI Marked Cube Mould (Cast Iron -for Concrete) 150mm Clamp Type (18 KG) Four-Part
- Cube Mould (Cast Iron - for Concrete) 150mm Clamp Type (16 KG) Two-Part

Apart from Cast Iron and Mild Steel Cube Moulds, we also manufacture Plastic & Polyurethane Cube Moulds



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265