

**Product Name :**  
Cylindrical Mould

**Product Code :**  
ASTConcrete-0005



**Description :**

Cylindrical Mould

**Technical Specification :**

The cylindrical mould is split vertically into two parts. The mean internal diameter is within  $\pm 0.2\text{mm}$  and height is within  $\pm 1\text{mm}$ . The ends are machined to  $\pm 0.05$ . The base plate and top plate are machined flat to  $\pm 0.03\text{mm}$ .

- Mould cylinder, cast iron, 150mm dia x 150mm height.
- Mould cylinder, cast iron, 150mm dia x 300mm height.
- Mould cylinder, cast iron, 160mm dia x 320mm height.
- Mould cylinder, cast iron, 100mm dia x 200mm height.
- Mould cylinder, cast iron, 300mm dia x 600mm height.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265