


Product Name :
Cement samplers

Product Code :
ASTCEMENT-0001



Description :

Cement samplers

Technical Specification :

Packaged cement tube sampler

Used to sample cement from packages. Made of brass, it has an outside dia. of 32 mm x 700 mm long.

- Weight approx. 2 kg

Bulk cement sampler

Used to sample cement in bulk storage or bulk shipments. It consists of two brass concentric tubes with slots. The inner tube rotates to close the slots and take the sample. The capacity of the internal tube is approx. 3 litres.

- Weight approx. 5 kg



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265