

Product Name :
Le Chatelier Moulds

Product Code :
ASTCEMENT-0003



Description :

Le Chatelier Moulds

Technical Specification :

Introduction

Used for determining the expansion of cement. The mould consists of a spring tensioned split cylinder 30 mm internal diameter, 30 mm high with two indicator stems which measure 165 mm from the points of the centre line of the cylinder and O ring. Two or three moulds are required for each test. To perform the test, a water bath is also required. See accessories.

We propose three packages:

Le Chatelier mould

Single unit. Individually checked, supplied complete with certificate of conformity and identified by a serial number.

- Weight approx.: 30 g

Le Chatelier mould. Pack of 6 units

Same as 62-L0025.

- Weight approx.: 150 g

Le Chatelier soundness kit

The kit includes all accessories to perform the test and to verify the conformity of the moulds.

Consisting of:

- Le Chatelier moulds
- Glass plates, 50 mm sq.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265