

**Product Name :**  
Sample Splitter

**Product Code :**  
ASTAGGRE-0001



**Description :**

Sample Splitter

**Technical Specification :**

Specification: Designed for the reduction of test samples which are too large in volume to be conveniently handled. It handles any material from Sand sizes up to dia. 108 mm. Each Chute bar is 12 mm wide so that openings of 12-24-36-48-60-72-84-96-108 mm are possible. Supplied complete with two collecting pans. Clam shell hopper: 30 litres capacity. Very sturdily constructed, it is totally Cadmium Plated for rust protection.

Applications    sampling and sample dividing

Field of application    agriculture, biology, chemistry / plastics, construction materials, engineering / electronics, environment / recycling, food, geology / metallurgy, glass / ceramics, medicine / pharmaceuticals

Feed material    bulk materials

Material feed size\*    RT 6.5: ~ 4 mm

RT 12.5: ~ 8 mm

RT 25: ~ 16 mm

RT 37.5: ~ 25 mm

RT 50: ~ 33 mm

RT 75: ~ 50 mm

Batch size / feed quantity\* RT 6.5: ? 3 l

RT 12.5: ? 3 l

RT 25: ? 16 l

RT 37.5: ? 16 l

RT 50: ? 16 l

RT 75: ? 16 l

Number of divisions 2

W x H x D RT 6.5: 30 x 27 x 25 cm

RT 12.5: 30 x 27 x 25 cm

RT 25: 62 x 42 x 26 cm

RT 37.5: 62 x 42 x 26 cm

RT 50: 62 x 42 x 26 cm

RT 75: 62 x 42 x 26 cm

Net weight RT 6.5: ~ 3.5 kg

RT 12.5: ~ 3.5 kg

RT 25: ~ 21.5 kg

RT 37.5: ~ 21.5 kg

RT 50: ~ 21.5 kg

RT 75: ~ 21.5 kg

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

