

Product Name :
Fill Box Apparatus (Kajima Test)

Product Code :
ASTConcrete-0023



Description :

Fill Box Apparatus (Kajima Test)

Technical Specification :

General descriptions: Fill Box Apparatus is used to measure the filling ability of self-compacting concrete with a maximum aggregate size of 20 mm. The apparatus is also known as “Kajima Test” Apparatus consists of a container (transparent) with a flat and smooth surface.

Dimensions 500x300x900 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265