


**Product Name :**  
Pocket Concrete Penetrometer

**Product Code :**  
ASTConcrete-0039



**Description :**

Pocket Concrete Penetrometer

**Technical Specification :**

Specification:

Consists of a needle having face area  $3/10$  sq. cm. and graduated at a distance of 25 cm. The needles point is an integral part of barrel which houses a calibrated spring. The spring is confined in a sleeve. The resistance offered by the concrete mortar is shown on the direct reading scale with a marker ring which holds its position when released. Scale range is 0-50 kg/cm<sup>2</sup> when the penetration resistance reaches a value of 35kg/cm<sup>2</sup> the concrete is assumed initially set. Supplied complete in carrying case.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

