

**Product Name :**  
Bleeding Of Freshly Mixed Concrete

**Product Code :**  
ASTConcrete-0040



**Description :**

Bleeding Of Freshly Mixed Concrete

**Technical Specification :**

Specification: This is used for determination of the relative quantity of mixing water that will bleed from a sample of freshly mixed concrete. Stainless steel cylindrical container for bleeding test to EN 480-4, approx. 14 liters capacity.

Internal dia 255 +/- 5mm x 280 +/- 5mm inside height. Complete with stainless steel cover. Flat rounded steel float with a diameter of (100 ± 10) mm.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

