

**Product Name :**  
RECIPROCATING PUMP DEMONSTRATION UNIT, Computer  
Interface

**Product Code :**  
FM177



**Description :**

A bench top unit with a reciprocating pump, 0.37kW motor with reduction gear, storage tank measuring tank, pressure relief valve, pulse chambers, sensors, and software. To be used with HP 100, HP 107 and computer or HP 104 and computer.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265