

**Product Name :**  
COMPACT RECIPROCATING PUMP TEST SET

**Product Code :**  
FM178



**Description :**

A bench top unit with reciprocating pump, 0.37kW motor, storage tank, pump pressure relief valve, pulse chambers, and measuring instruments. Optional : Advanced inverter for speed control as well as for motor speed and output power indication.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265