

Product Name :
COMPACT RECIPROCATING PUMP TEST SET

Product Code :
FM178



Description :

A bench top unit with reciprocating pump, 0.37kW motor, storage tank, pump pressure relief valve, pulse chambers, and measuring instruments. Optional : Advanced inverter for speed control as well as for motor speed and output power indication.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265