

**Product Name :**  
Load Ring Penetrometer

**Product Code :**  
ASTSOIL-0028



**Description :**

Load Ring Penetrometer

**Technical Specification :**

Specification: Used for measuring the bearing strength and compaction degree of soils. The apparatus consists of a "T" shaped handle connected to a load ring 1 kN (100 kgf) cap., with max load pointer, and of an extension rod with five 100 mm graduations. The 30° end cone has an area of 645 mm<sup>2</sup> (1 sq. in). Supplied complete with calibration chart.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265