


Product Name :
Load Ring Penetrometer

Product Code :
ASTSOIL-0028



Description :

Load Ring Penetrometer

Technical Specification :

Specification: Used for measuring the bearing strength and compaction degree of soils. The apparatus consists of a "T" shaped handle connected to a load ring 1 kN (100 kgf) cap., with max load pointer, and of an extension rod with five 100 mm graduations. The 30° end cone has an area of 645 mm² (1 sq. in). Supplied complete with calibration chart.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265