


Product Name :
Pocket Penetrometers

Product Code :
ASTSOIL-0030



Description :

Pocket Penetrometers

Technical Specification :

Specification:

Designed for making field classification of cohesive soils in terms of consistency, shear strength and approximate unconfined compressive strength.

- Stainless steel construction.
- Tip Diameters: 4.5 mm dia. for very hard soil, 6.35 mm for medium and soft soil, 8.98 mm for soft soil.
- Measuring range: 0 to 1000 kPa
- Dimensions (assembled): 210 mm length x 20 mm dia. approx.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265